

# Slim Handpiece

Appropriate for MIS  
biportal endoscopic spine surgeries




Primado Series

# Slim Handpiece 200

## Enhanced Accessibility and Durability

Features sleeve bearings tailored to meet the demands of endoscopic spine surgeries in aqueous environments

- Improved visibility and operability at the surgical site with a reduced tip diameter
- A simplified design minimizes the risk of post-surgical mechanical failure



Description	Order Code	Dimension of Bur Guard (mm)	Maximum Rotation Speed (min <sup>-1</sup> )
Slim Handpiece200 155/U	P200-1A155-U	Φ2.1-3.3×L100.8	80,000

\*Included product: Wrench for Slim/Super Slim (P200-1ARA-W)

## Applicable Motors/Cutting Accessories



Order Code  
**P200-SMH-S**  
Primado2 Slim Motor Short



Order Code  
**P200-SMH**  
Primado2 Slim Motor



Order Code  
**P200-SMH-HS**  
Primado2 Slim Motor/HS  
• Hand Control Type



Order Code  
**PDA-SMH**  
Primado Air Slim Motor



Order Code  
**PDA-SMH-HS**  
Primado Air Slim Motor/HS  
• Hand Control Type



**Diamond Bur**  
Head Size (mm) **2.0/3.0/4.0**



**Coarse Diamond Bur**  
Head Size (mm) **2.0/3.0/4.0**



**Neuro Diamond Bur**  
Head Size (mm) **2.0/3.0/4.0**



**Neuro Coarse Diamond Bur**  
Head Size (mm) **2.0/3.0/4.0**



**Round Fluted Bur**  
Head Size (mm) **2.0/3.0/4.0**



**Match Head Cutting Bur**  
Head Size (mm) **2.5/3.0**